

## Product Datasheet

# Fab Golden Syrian Hamster IgG (H&L) Antibody Rhodamine Conjugated (orb348445)

|                            |                                                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Fab Golden Syrian Hamster IgG (H&L) antibody (TRITC)                                                                                                                                                                                                                                                                                                                                              |
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Reactivity</b>          | Golden Hamster                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Conjugation</b>         | TRITC                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Tested Applications</b> | FC, FLISA, IF                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen</b>           | Golden Syrian Hamster IgG whole molecule                                                                                                                                                                                                                                                                                                                                                          |
| <b>Preservatives</b>       | 0.01% (w/v) Sodium Azide                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Concentration</b>       | 1.0 mg/mL                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>             | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Application notes</b>   | This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. |
| <b>Isotype</b>             | IgG Fab                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                        |

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antibody Type</b>   | Secondary Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Purity</b>          | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Golden Syrian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein and anti-Goat Serum. No reaction was observed against anti-Papain or anti-Goat IgG F(c). |
| <b>Dilution Range</b>  | FLISA: 1:10,000 - 1:50,000, FC: 1:500 - 1:2,500, IF: 1:1,000 - 1:5,000                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                                                  |

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)