

Product Datasheet

Major allergen Api g 1, isoallergen 2 Antibody, HRP conjugated (orb350933)

Catalog Number	orb350933
Category	Antibodies
Description	This Major allergen Api g 1, isoallergen 2 Antibody, HRP conjugated is an HRP polyclonal product. It targets Major allergen Api g 1, isoallergen 2 using a Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified.
Target	Major allergen Api g 1, isoallergen 2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	HRP
Reactivity	Plant
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA)
UniProt ID	P92918
Tested applications	ELISA

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)