

## Product Datasheet

### Major allergen Api g 1, isoallergen 2 Antibody, FITC conjugated (orb350934)

<b>Catalog Number</b>	orb350934
<b>Category</b>	Antibodies
<b>Description</b>	This Major allergen Api g 1, isoallergen 2 Antibody, FITC conjugated is an FITC polyclonal product. It targets Major allergen Api g 1, isoallergen 2 using a Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA) as the immunogen. Purification: >95%, Protein G purified.
<b>Target</b>	Major allergen Api g 1, isoallergen 2
<b>Clonality</b>	Polyclonal
<b>Species/Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Conjugation</b>	FITC
<b>Reactivity</b>	Plant
<b>Form/Appearance</b>	Liquid
<b>Buffer/Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Purification</b>	>95%, Protein G purified
<b>Immunogen</b>	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA)
<b>UniProt ID</b>	<b>P92918</b>
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)