

Product Datasheet

GSTK1 Antibody (HRP) (orb352214)

| | |
|-----------------------------|---|
| Catalog Number | orb352214 |
| Category | Antibodies |
| Description | This GSTK1 Antibody, HRP conjugated is an HRP polyclonal product. It targets GSTK1 using a Recombinant Human Glutathione S-transferase kappa 1 protein (2-226AA) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified. |
| Target | GSTK1 |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | HRP |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Buffer/Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Purification | >95%, Protein G purified |
| Immunogen | Recombinant Human Glutathione S-transferase kappa 1 protein (2-226AA) |
| UniProt ID | Q9Y2Q3 |
| Tested applications | ELISA |
| Dilution range | ELISA |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)