

Product Datasheet

ACTRT1 Antibody, HRP conjugated (orb353806)

Catalog Number	orb353806
Category	Antibodies
Description	This ACTRT1 Antibody, HRP conjugated is an HRP polyclonal product. It targets Major allergen Api g 1, isoallergen 2 using a Recombinant Human Actin-related protein T1 protein (201-376AA) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified.
Target	Major allergen Api g 1, isoallergen 2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	HRP
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Actin-related protein T1 protein (201-376AA)
UniProt ID	Q8TDG2
Tested applications	ELISA
Dilution range	WB:1:500-5000, IF:1:50-200

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)