

Product Datasheet

Rabbit Goat IgG Antibody, Biotin conjugated (orb357962)

| | |
|-----------------------------|--|
| Catalog Number | orb357962 |
| Category | Antibodies |
| Description | This Rabbit anti-Goat IgG Antibody;Biotin conjugated is an Biotin polyclonal product, using a IgG purified from goat serum as the immunogen. This antibody is suitable for ELISA. Purification: >95%, ProteinG purified. |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | Biotin |
| Reactivity | Goat |
| Form/Appearance | Liquid |
| Buffer/Preservatives | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification | >95%, ProteinG purified |
| Immunogen | IgG purified from goat serum |
| Tested applications | ELISA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)