

Product Datasheet

Rat Q91Y39 protein (Active) (orb359141)

Catalog Number	orb359141
Category	Proteins
Description	This Rat Q91Y39 protein (Active) spans the amino acid sequence from region 20-135aa. Purity: > 96% as determined by SDS-PAGE and HPLC.
Tag	Tag-Free
Form/Appearance	Lyophilized powder
Buffer/Preservatives	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Purity	> 96% as determined by SDS-PAGE and HPLC.
Protein Sequence	SEAILPIASSCCTEVSHHIPRLLERVNSCSIQRADGDCDLAAVILHVKRRRICVSPHNPTLKR WMSASEMKNKGENLCPRKKQDSGKDRKGHTPRKHGKHGTRRIHGTHDHEAPR
Protein Length	Full Length of Mature Protein
UniProt ID	Q91Y39
MW	13.1 kDa
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Endotoxins	Less than 1.0 EU/µg as determined by LAL method.
Source	E.Coli
Biological Origin	Rattus norvegicus (Rat)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biological Activity

Fully biologically active when compared to standard. The biologically active determined by a chemotaxis bioassay using human lymphocytes is in a concentration range of 5.0-50 ng/ml.

Expression Region

20-135aa

Storage


The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Note

For research use only

Expiration Date

6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)