

Product Datasheet

NR1I2 Rabbit Polyclonal Antibody (orb385455)

Catalog Number	orb385455
Category	Antibodies
Description	NR1I2 Rabbit Polyclonal Antibody
Target	NR1I2. This antibody has specificity with NR1I2.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Porcine, Rat
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Immunogen	KLH conjugated synthetic peptide derived from human NR1I2. Please contact us for the exact immunogen sequence. The peptide is available as orb13108.
UniProt ID	O54915, O75469, Q9R1A7
MW	50 kDa
Tested applications	ICC, IF, IHC-P, WB
Dilution range	IHC-P: 1:100-600, IF/ICC: 1:100-600, WB: 1:200-1500

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Entrez	8856
NCBI	003889, 003880
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin embedded pig kidney tissue using anti-NR1I2 (2.5 ug/ml)



Immunohistochemical staining of paraffin embedded pig kidney tissue using NR1I2 antibody (2.5 ug/ml)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



IHC-P image of human breast cancer tissue using NR112 antibody (2.5 ug/ml)



WB analysis of mouse liver (lane 1), rat kidney (lane 2), rat colon (lane 3), rat small intestines (lane 4), human breast cancer (lane 5); hela cells (lane 6) using anti-NR112 (1 ug/ml)



IHC-P staining of rat colon tissue using anti-NR112 (2.5 ug/ml)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)