

Product Datasheet

Cyclin D1 Antibody / CCND1 (orb386261)


Description	Mouse monoclonal antibody against Cyclin D1 conjugated to PE
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FACS
Immunogen	Human recombinant full length Cyclin D1 protein was used as the immunogen.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Cyclin D1 antibody at 2-8°C. Conjugate is light sensitive, store in the dark.
Note	For research use only
Application notes	The dilution stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Cyclin D1 antibody to be titered up or down for optimal performance.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Clone Number	BCLD1-1
Antibody Type	Primary Antibody
Purity	Protein G affinity chromatography
Hazard Information	This Cyclin D1 antibody is available for research use only.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Entrez **595****Dilution Range** Flow Cytometry: 5ul/test/million cells or 5ul/test/100ul of whole blood**Expiration Date** 12 months from date of receipt.

FACS testing of human HeLa cells: Black = cells alone; Green = isotype control; Red = Cyclin D1 antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)