

Product Datasheet

Beta-Galactosidase Activity Colorimetric Microplate Assay Kit (orb390805)

Catalog Number	orb390805
Category	Assays and Kits
Description	Detection and Quantification of Beta-Galactosidase Activity.
Reactivity	All
Tested applications	FA
Application notes	A spectrophotometer can be used, but a micro cuvette with a volume of $\leq 200 \mu\text{L}$ is required.
Assay Type	Detection
Range	10 $\mu\text{mol/L}$ - 1000 $\mu\text{mol/L}$
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9031-11-2
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Representative Standard Curve



Positive Control reaction in 96-well plate assay with decreasing the concentration

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558