


Product Datasheet

Serum Kalium Microplate Assay Kit (orb390810)

Catalog Number	orb390810
Category	Assays and Kits
Description	Detection and Quantification of Serum Kalium Content.
Reactivity	All
Tested applications	FA
Assay Type	Detection
Range	0.2 mmol/L - 2 mmol/L
Sample Types	Serum, Urine, Saliva, Other biological fluids
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Expiration Date	12 months from date of receipt.



Representative Standard Curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558