

## Product Datasheet

# PRKCI (Phospho-T412/410) Blocking Peptide (orb392840)

<b>Catalog Number</b>	orb392840
<b>Category</b>	Proteins
<b>Description</b>	The peptide is used to block Anti-PKC iota/zeta (pT412/410) Antibody (orb393138) reactivity.
<b>Target</b>	PRKCI; PRKCZ
<b>Reactivity</b>	Canine, Human, Mouse, Rabbit, Rat, Zebrafish
<b>Form/Appearance</b>	Lyophilized powder
<b>Purity</b>	>85%
<b>UniProt ID</b>	<b>Q02956, P09217, Q62074, P41743, Q05513, F1M7Y5</b>
<b>Tested applications</b>	Blocking
<b>Application notes</b>	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Source</b>	Synthetic
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Entrez</b>	<b>5590, 18762, 25522, 18759, 5584, 84006</b>
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)