

## Product Datasheet

### HLA-C Antibody, FITC conjugated (orb401646)

|                        |                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------|
| <b>Description</b>     | HLA-C Antibody, FITC conjugated                                                                |
| <b>Species/Host</b>    | Rabbit                                                                                         |
| <b>Reactivity</b>      | Human                                                                                          |
| <b>Conjugation</b>     | FITC                                                                                           |
| <b>Immunogen</b>       | Recombinant Human HLA class I histocompatibility antigen, Cw-3 alpha chain protein (200-306AA) |
| <b>Target</b>          | HLA-C                                                                                          |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                                    |
| <b>Form/Appearance</b> | Liquid                                                                                         |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                  |
| <b>Note</b>            | For research use only                                                                          |
| <b>Isotype</b>         | IgG                                                                                            |
| <b>Clonality</b>       | Polyclonal                                                                                     |
| <b>Uniprot ID</b>      | <b>P04222</b>                                                                                  |
| <b>Dilution Range</b>  | ELISA                                                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)