

Product Datasheet

Chicken IFN gamma ELISA Kit (orb403361)

Catalog Number	orb403361
Description	Chicken Interferon gamma; ,IFN-gamma; ELISA Kit
Reactivity	Gallus
Range	3.125 pg/mL-800 pg/mL
Target	IFNG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Assay Time	1-5h
Uniprot ID	P49708
Sensitivity	3.125 pg/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)