

## Product Datasheet

### Chicken growth hormone ELISA Kit (orb403416)

<b>Catalog Number</b>	orb403416
<b>Category</b>	Assays and Kits
<b>Description</b>	The Chicken growth hormone ELISA Kit detects the target GH in gallus samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 625 pg/mL-10000 pg/mL. The kit offers a sensitivity of 625 pg/mL.
<b>Target</b>	GH
<b>Reactivity</b>	Gallus
<b>UniProt ID</b>	<b>P08998</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Competitive
<b>Assay Time</b>	1-5h
<b>Range</b>	625 pg/mL-10000 pg/mL
<b>Sensitivity</b>	625 pg/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)