

Product Datasheet

Dog TSH ELISA Kit (orb403534)

Catalog Number	orb403534
Category	Assays and Kits
Description	The Dog TSH ELISA Kit detects the target TSH in canine samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.625 μ U/mL-10 μ U/mL. The kit offers a sensitivity of 0.625 μ U/mL.
Target	TSH
Reactivity	Canine
UniProt ID	P54828
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.625 μ U/mL-10 μ U/mL
Sensitivity	0.625 μ U/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)