

Product Datasheet

Goat CGA ELISA Kit (orb403695)

Catalog Number	orb403695
Category	Assays and Kits
Description	The Goat CGA ELISA Kit detects the target CGA in goat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.5 μ U/mL-20 μ U/mL. The kit offers a sensitivity of 0.25 μ U/mL.
Target	CGA
Reactivity	Goat
UniProt ID	Q8WMW8
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.5 μ U/mL-20 μ U/mL
Sensitivity	0.25 μ U/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)