

## Product Datasheet

### Goat Luteinizing Hormone ELISA Kit (orb403709)

<b>Catalog Number</b>	orb403709
<b>Category</b>	Assays and Kits
<b>Description</b>	The Goat Luteinizing Hormone ELISA Kit detects the target LH in goat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.3 mLU/mL-60 mLU/mL. The kit offers a sensitivity of 0.15 mIU/mL.
<b>Target</b>	LH
<b>Reactivity</b>	Goat
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Competitive
<b>Assay Time</b>	1-5h
<b>Range</b>	0.3 mLU/mL-60 mLU/mL
<b>Sensitivity</b>	0.15 mIU/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)