

Product Datasheet

Goat TNF ELISA Kit (orb403732)

Catalog Number	orb403732
Category	Assays and Kits
Description	The Goat TNF ELISA Kit detects the target TNF in goat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 pg/mL-1000 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
Target	TNF
Reactivity	Goat
UniProt ID	P13296
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)