

Product Datasheet

Guinea Pig OVA sIgE ELISA Kit (orb403774)

Catalog Number	orb403774
Category	Assays and Kits
Description	The Guinea Pig OVA sIgE ELISA Kit detects the target OVA sIgE in guinea pig samples. Sample types including serum and plasma.
Target	OVA sIgE
Reactivity	Guinea pig
Application notes	Detection Wavelength: 450 nm
Assay Type	Indirect
Assay Time	1-5h
Range	Request Information
Sensitivity	Request Information
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)