

## Product Datasheet

### Guinea Pig IgE ELISA Kit (orb403777)

<b>Catalog Number</b>	orb403777
<b>Category</b>	Assays and Kits
<b>Description</b>	The Guinea Pig IgE ELISA Kit detects the target IgE in guinea pig samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 93.75 ng/mL-6000 ng/mL. The kit offers a sensitivity of 23.4 ng/mL.
<b>Target</b>	IgE
<b>Reactivity</b>	Guinea pig
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	93.75 ng/mL-6000 ng/mL
<b>Sensitivity</b>	23.4 ng/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)