

Product Datasheet

Human ARG1 ELISA Kit (orb403993)

Catalog Number	orb403993
Category	Assays and Kits
Description	The Human ARG1 ELISA Kit detects the target ARG1 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and ascitic fluid. The detection range is 1.56 ng/mL-100 ng/mL. The kit offers a sensitivity of 0.39 ng/mL.
Target	ARG1
Reactivity	Human
UniProt ID	P05089
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Sample Types	serum, plasma, tissue homogenates, ascitic fluid
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)