

Product Datasheet

Human CSPG4 ELISA Kit (orb404145)

Catalog Number	orb404145
Category	Assays and Kits
Description	The Human CSPG4 ELISA Kit detects the target CSPG4 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 25 pg/mL-1600 pg/mL. The kit offers a sensitivity of 6.25 pg/mL.
Target	CSPG4
Reactivity	Human
UniProt ID	Q6UVK1
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)