

Product Datasheet

Human GRIN1 ELISA Kit (orb404487)

Catalog Number	orb404487
Category	Assays and Kits
Description	The Human GRIN1 ELISA Kit detects the target GRIN1 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 125 pg/mL-8000 pg/mL. The kit offers a sensitivity of 31.25 pg/mL.
Target	GRIN1
Reactivity	Human
UniProt ID	Q05586
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)