

Product Datasheet

Human IL22RA2 ELISA Kit (orb404664)

Catalog Number	orb404664
Category	Assays and Kits
Description	The Human IL22RA2 ELISA Kit detects the target IL22RA2 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.81 pg/mL.
Target	IL22RA2
Reactivity	Human
UniProt ID	Q969J5
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)