

Product Datasheet

Human LSR ELISA Kit (orb404750)

Catalog Number	orb404750
Category	Assays and Kits
Description	The Human LSR ELISA Kit detects the target LSR in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 6.25 pg/mL-400 pg/mL. The kit offers a sensitivity of 1.56 pg/mL.
Target	LSR
Reactivity	Human
UniProt ID	Q86X29
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)