

## Product Datasheet

### Human MSN ELISA Kit (orb404850)

<b>Catalog Number</b>	orb404850
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human MSN ELISA Kit detects the target MSN in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.81 pg/mL.
<b>Target</b>	MSN
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>P26038</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Sensitivity</b>	7.81 pg/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)