

Product Datasheet

Human SERPINB2 ELISA Kit (orb405027)

Catalog Number	orb405027
Category	Assays and Kits
Description	The Human SERPINB2 ELISA Kit detects the target SERPINB2 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates and cell lysates. The detection range is 15.6 pg/mL-1000 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
Target	SERPINB2
Reactivity	Human
UniProt ID	P05120
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)