

Product Datasheet

Human PDPN ELISA Kit (orb405035)

Catalog Number	orb405035
Category	Assays and Kits
Description	The Human PDPN ELISA Kit detects the target PDPN in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 25 pg/mL-1600 pg/mL. The kit offers a sensitivity of 6.25 pg/mL.
Target	PDPN
Reactivity	Human
UniProt ID	Q86YL7
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)