

## Product Datasheet

### Human PDCD1LG2 ELISA Kit (orb405069)

<b>Catalog Number</b>	orb405069
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human PDCD1LG2 ELISA Kit detects the target PDCD1LG2 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.45 ng/ml-30 ng/ml. The kit offers a sensitivity of 0.11 ng/ml.
<b>Target</b>	PDCD1LG2
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>Q9BQ51</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	0.45 ng/ml-30 ng/ml.
<b>Sensitivity</b>	0.11 ng/ml
<b>Sample Types</b>	serum, plasma, tissue homogenates.
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)