

Product Datasheet

Human RTN4 ELISA Kit (orb405226)

| | |
|--------------------------|---|
| Catalog Number | orb405226 |
| Category | Assays and Kits |
| Description | The Human RTN4 ELISA Kit detects the target RTN4 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 46.88 pg/mL-3000 pg/mL. The kit offers a sensitivity of 11.75 pg/mL. |
| Target | RTN4 |
| Reactivity | Human |
| UniProt ID | Q9NQC3 |
| Application notes | Detection Wavelength: 450 nm |
| Assay Type | Sandwich |
| Assay Time | 1-5h |
| Range | 46.88 pg/mL-3000 pg/mL |
| Sensitivity | 11.75 pg/mL |
| Sample Types | serum, plasma, tissue homogenates |
| Note | For research use only |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)