

## Product Datasheet

### Human WISP2 ELISA Kit (orb405555)

<b>Catalog Number</b>	orb405555
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human WISP2 ELISA Kit detects the target WISP2 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.156 ng/mL-10 ng/mL. The kit offers a sensitivity of 0.039 ng/mL.
<b>Target</b>	WISP2
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>O76076</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Sensitivity</b>	0.039 ng/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)