

Product Datasheet

Human YWHAE ELISA Kit (orb405581)

Catalog Number	orb405581
Category	Assays and Kits
Description	The Human YWHAE ELISA Kit detects the target YWHAE in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 78 pg/mL-5000 pg/mL. The kit offers a sensitivity of 19.5 pg/mL.
Target	YWHAE
Reactivity	Human
UniProt ID	P62258
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	78 pg/mL-5000 pg/mL
Sensitivity	19.5 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)