

Product Datasheet

Human CD109 ELISA Kit (orb405670)

Catalog Number	orb405670
Category	Assays and Kits
Description	The Human CD109 ELISA Kit detects the target CD109 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 94 pg/mL-6000 pg/mL. The kit offers a sensitivity of 23.4 pg/mL.
Target	CD109
Reactivity	Human
UniProt ID	Q6YHK3
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	94 pg/mL-6000 pg/mL
Sensitivity	23.4 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)