

Product Datasheet

Human KRAS ELISA Kit (orb405742)

Catalog Number	orb405742
Category	Assays and Kits
Description	The Human KRAS ELISA Kit detects the target KRAS in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 23.5 pg/mL-1500 pg/mL. The kit offers a sensitivity of 5.8 pg/mL.
Target	KRAS
Reactivity	Human
UniProt ID	P01116
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)