

Product Datasheet

Human NMDA antibody (IgG) ELISA Kit (orb405773)

Catalog Number	orb405773
Description	Human N-methyl-D-aspartic acid (NMDA) antibody (IgG) ELISA Kit
Reactivity	Human
Range	Request Information
Target	NMDA Ab (IgG)
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum
Assay Time	1-5h
Sensitivity	Request Information
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)