

Product Datasheet

Human ENO1 ELISA Kit (orb405903)

Catalog Number	orb405903
Category	Assays and Kits
Description	The Human ENO1 ELISA Kit detects the target ENO1 in human samples. It measures levels in sample types including serum and tissue homogenates and plasma. The detection range is 0.312 ng/mL-20 ng/mL. The kit offers a sensitivity of 0.078 ng/mL.
Target	ENO1
Reactivity	Human
UniProt ID	P06733
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Sample Types	serum, tissue homogenates, plasma
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)