

Product Datasheet

Human amyloid beta peptide 1-40 ELISA Kit (orb405915)

Catalog Number	orb405915
Category	Assays and Kits
Description	The Human amyloid beta peptide 1-40 ELISA Kit detects the target A β 1-40 in human samples. It measures levels in sample types including serum and tissue homogenates and cerebrospinal fluid (CSF). The detection range is 125 pg/mL-8000 pg/mL. The kit offers a sensitivity of 31.25 pg/mL.
Target	A β 1-40
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Sample Types	serum, tissue homogenates, cerebrospinal fluid (CSF)
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)