

## Product Datasheet

### Human IL17C ELISA Kit (orb405972)

|                          |  |
|--------------------------|--|
| <b>Catalog Number</b>    | orb405972  |
| <b>Category</b>          | Assays and Kits  |
| <b>Description</b>       | The Human IL17C ELISA Kit detects the target IL17C in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.156 ng/ml -10 ng/ml. The kit offers a sensitivity of 0.039 ng/ml. |
| <b>Target</b>            | IL17C  |
| <b>Reactivity</b>        | Human  |
| <b>UniProt ID</b>        | <b>Q9P0M4</b>  |
| <b>Application notes</b> | Detection Wavelength: 450 nm   |
| <b>Assay Type</b>        | Sandwich   |
| <b>Assay Time</b>        | 1-5h   |
| <b>Range</b>             | 0.156 ng/ml -10 ng/ml  |
| <b>Sensitivity</b>       | 0.039 ng/ml  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates  |
| <b>Note</b>              | For research use only  |
| <b>Expiration Date</b>   | Please enquire.  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)