

Product Datasheet

Human IL18BP ELISA Kit (orb405975)

Catalog Number	orb405975
Category	Assays and Kits
Description	The Human IL18BP ELISA Kit detects the target IL18BP in human samples. It measures levels in sample types including serum and plasma. The detection range is 15.6 pg/ml-1000 pg/ml. The kit offers a sensitivity of 3.9 pg/ml.
Target	IL18BP
Reactivity	Human
UniProt ID	O95998
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)