

Product Datasheet

Human IL1 Receptor ELISA Kit (orb405981)

Catalog Number	orb405981
Category	Assays and Kits
Description	The Human IL1 Receptor ELISA Kit detects the target IL1R1 in human samples. It measures levels in sample types including serum and plasma and tissue homogenates. The detection range is 0.156 ng/ml-10 ng/ml. The kit offers a sensitivity of 0.039 ng/ml.
Target	IL1R1
Reactivity	Human
UniProt ID	P14778
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Sample Types	serum, plasma,tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)