

Product Datasheet

Human IL-28B ELISA Kit (orb406001)

Catalog Number	orb406001
Category	Assays and Kits
Description	The Human IL-28B ELISA Kit detects the target IL-28B in human samples. It measures levels in sample types including serum and plasma and cell culture supernats and tissue homogenates. The detection range is 1.56 pg/mL-100 pg/mL. The kit offers a sensitivity of 0.39 pg/mL.
Target	IL-28B
Reactivity	Human
UniProt ID	Q8IZI9
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.56 pg/mL-100 pg/mL
Sensitivity	0.39 pg/mL
Sample Types	serum, plasma, cell culture supernats, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)