

Product Datasheet

Human PON1 ELISA Kit (orb406214)

Catalog Number	orb406214
Category	Assays and Kits
Description	The Human PON1 ELISA Kit detects the target PON1 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 31.25 mIU/mL-2000 mIU/mL. The kit offers a sensitivity of 7.81 mIU/mL.
Target	PON1
Reactivity	Human
UniProt ID	P27169
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 mIU/mL-2000 mIU/mL
Sensitivity	7.81 mIU/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)