

Product Datasheet

Human CNN2 ELISA Kit (orb406286)

Catalog Number	orb406286
Category	Assays and Kits
Description	The Human CNN2 ELISA Kit detects the target CNN2 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.8 pg/mL.
Target	CNN2
Reactivity	Human
UniProt ID	Q99439
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)