

## Product Datasheet

### Human BMP15 ELISA Kit (orb406350)

<b>Catalog Number</b>	orb406350
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human BMP15 ELISA Kit detects the target BMP15 in human samples. It measures levels in sample types including serum, plasma and cell culture supernates. The detection range is 15.6 pg/mL-1000 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
<b>Target</b>	BMP15
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>O95972</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Sensitivity</b>	3.9 pg/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)