

## Product Datasheet

### Human PRDX2 ELISA Kit (orb406390)

<b>Catalog Number</b>	orb406390
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human PRDX2 ELISA Kit detects the target PRDX2 in human samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 47 pg/mL-3000 pg/mL. The kit offers a sensitivity of 11.75 pg/mL.
<b>Target</b>	PRDX2
<b>Reactivity</b>	Human
<b>UniProt ID</b>	<b>P32119</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	47 pg/mL-3000 pg/mL
<b>Sensitivity</b>	11.75 pg/mL
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)