

Product Datasheet

Human ALOX12 ELISA Kit (orb406428)

Catalog Number	orb406428
Category	Assays and Kits
Description	The Human ALOX12 ELISA Kit detects the target ALOX12 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 1.25 ng/mL-80 ng/mL. The kit offers a sensitivity of 0.31 ng/mL.
Target	ALOX12
Reactivity	Human
UniProt ID	P18054
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.25 ng/mL-80 ng/mL
Sensitivity	0.31 ng/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)