

Product Datasheet

Human SPAG9 ELISA Kit (orb406541)

Catalog Number	orb406541
Category	Assays and Kits
Description	The Human SPAG9 ELISA Kit detects the target SPAG9 Ab in human samples. It measures levels in sample types including serum and plasma. The detection range is 50 mU/mL-1200 mU/mL. The kit offers a sensitivity of 20 mU/mL.
Target	SPAG9 Ab
Reactivity	Human
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	50 mU/mL-1200 mU/mL
Sensitivity	20 mU/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)