

Product Datasheet

Human Galectin 3 ELISA Kit (orb406712)

Catalog Number	orb406712
Category	Assays and Kits
Description	The Human Galectin 3 ELISA Kit detects the target LGALS3 in human samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates and saliva. The detection range is 0.625 ng/mL-40 ng/mL. The kit offers a sensitivity of 0.156 ng/mL.
Target	LGALS3
Reactivity	Human
UniProt ID	P17931
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, saliva
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)