

Product Datasheet

Human TNFRSF11A ELISA Kit (orb406716)

Catalog Number	orb406716
Category	Assays and Kits
Description	The Human TNFRSF11A ELISA Kit detects the target TNFRSF11A in human samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.312 ng/ml-20 ng/ml. The kit offers a sensitivity of 0.078 ng/ml.
Target	TNFRSF11A
Reactivity	Human
UniProt ID	Q9Y6Q6
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.3125 ng/ml- 20ng/ml
Sensitivity	0.078 ng/ml
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)